

Professional Devices

ComboCare™ E-Stim and Ultrasound Combo Professional Device

The ComboCare™ is the latest addition to our professional line of therapy systems. This all-new microprocessor controlled unit combines clinical ultrasound therapy with five treatment waveforms: IF 4-Pole, IF 2-Pole, Russian, EMS, and TENS. Large easy-to-read LCD panel allows the clinician to quickly set up treatment parameters. A dual frequency 5 cm sound head comes standard with the unit.



Product Features:

- Two channels of electrical stimulation output
- Large easy-to-read graphical LCD display with exact dial performance and programmable timer
- Waveforms include TENS, EMS, 2-Pole Interferential, 4-Pole Interferential and Russian
- Dual clinical ultrasound features a 5 cm sound head which offers 1 MHz and 3 MHz frequencies
- Ultrasound has 10% - 100% duty cycle selections
- 10 presets per waveform totaling 50 in all
- 20 ultrasound presets
- User defined memory positions for user protocols
- Designed for tabletop, therapy cart or mobile use
- Includes eight sets of electrodes, eight sets of rubber electrodes with sponges, eight straps, four sets of lead wires and a user manual
- Product dimensions: 10 ½" (W) x 15" (L) x 8" (H)
- Product weight: 6.7 lbs
- One year warranty

Product Specifications:

General Unit Specifications

Wavelength:	IFC-4P, IFC-2P, TENS, EMS and Russian
Pulse Width:	30-40 μ s
Frequency:	1 MHz and 3 MHz for ultrasound, 1-250 Hz for LF, beat frequency for IF: 1-150 Hz, carrier frequency for IF: 4KHz, carrier frequency for Russian: 2.5KHz
Power Supply:	Adaptor, input 100-240V 50/60 Hz, output 15V/3A
Intensity:	Max 3W/CM ² (pulse) for ultrasound, 0-50mA for IF (1K Ω), 0-100mA for LF (1K Ω)
Working Time:	Max 60 minutes for IF, Max 30 minutes for ultrasound and combo
Dimensions:	9.8" x 7.3" x 3.2"
Weight:	5.5 lbs.
Working Environment:	Temperature: 10°C - 40°C Humidity: 30% RH to 85% RH Atmospheric Pressure: 800hPa - 1060hPa
Storage Environment:	Temperature: 10°C - 55°C Humidity: 20% RH to 85% RH Atmospheric Pressure: 700hPa - 1060hPa

Item #	Description	Price
18-0450	ComboCare™ E-Stim and Ultrasound Combo Professional Device	\$1,295
32-0511	1 x 5 Litre Wavelength Multi-Purpose BLUE Ultrasound Gel	\$10.50
19-1090	BodyMed 2" x 2" Square Self-Adhering Electrodes (pack of 4 pcs)	\$4.95
19-1093	BodyMed 3" Round Self-Adhering Electrodes (pack of 4 pcs)	\$7.95

* Prescription Required: This is a Class II Medical Device. Health Canada restricts these devices to sale by or on the order of a licensed health care professional.

Kosim
ideal for health
www.kosim.ca

**To Order: (416) 502-3700, Toll Free: 1-877-KOSIM50
Fax: (416) 502-2808 or sales@kosim.ca**