


# Intelect<sup>®</sup>

## LEGEND<sup>XT</sup>




## An Exceptional Therapy System!

- Modular Design
- Patient Data Cards
  - Pain Profiles
  - User Defined Protocols



ELECTRODE PLACEMENT



CLINICAL WAVEFORMS



PAIN MAP PROFILES



USER PROTOCOLS


# An Exceptional Therapy System!

**Intelect is a name** the clinical rehabilitation community has relied on for many years. Well known for its exceptional clinical performance, continual quality, low cost and ease of use, Intelect has evolved to an exceptional new design and function. The new Intelect Legend XT is a clean modular design available in two or four channel electrotherapy or combination electrotherapy ultrasound systems. A modular design means affordability as it allows the clinician to select a two-channel system now and add two more channels later without significant cost.

**Clinical performance** has increased by providing six essential electrotherapy waveforms. Complete session parameters can be stored on Patient Data Cards. The new on-board Patient Data Management (PDM) system, along with the Patient Data Card, work together to overcome the problems of time management, the inconsistency of patient care and the increased need for accurate patient documentation.

**Comprehensive patient pain profiles** are recorded with a VAS (Visual Analog Scale), Numeric Pain Index and Pain Map Profile for each patient. Electrode placement diagrams can be created for each patient to improve the consistency of care. Session notes are also accessible for each therapy session. User defined protocols are another great feature to customize any waveform to your specific needs and store it by name.

**Attention to detail** is logically designed in the user interface and continues into the new therapy system cart, providing an integrated cord management strategy to minimize the tangling of cords. Six concealed storage bins keep therapy resources such as gels, towels and electrodes close at hand.



**Easy to Add Modularity** 2 more channels or a battery module

## PATIENT DOCUMENTATION

Allows clinicians to capture and record accurate patient information through the following programs: Pain Map Profile, Modality Records, Electrode Placement Diagrams and Session Notes.



## MULTIPLE WAVEFORMS

The essential clinical waveforms to manage a broad spectrum of clinical conditions. Interferential, Premodulated, High-Volt, Russian, Symmetrical Biphasic and Microcurrent.



## USER DEFINED PROTOCOLS

Customizes any waveform to your specific needs. User protocols are organized alphabetically by name for easy recall.



## FOUR ERGONOMIC APPLICATORS

Four different sizes for all clinician needs; fits naturally in the hand with an indicator light; Available in 2cm<sup>2</sup>, 5cm<sup>2</sup> and 10 cm<sup>2</sup> applicators; NEW 1 cm<sup>2</sup> applicator for hand therapy.



## PORTABILITY MADE EASY

The optional Battery Module is designed for treatment away from outlets, allowing therapy to be performed anywhere throughout the clinic. A specially-designed optional Carry Case can also be used to transport the Intelect Legend XT and its accessories.



# Intelect<sup>®</sup> LEGEND XT

For more information, downloads and an online demo, go to our website:

[www.intelect-legendxt.com](http://www.intelect-legendxt.com)

## DISPLAY

5" FSTN LCD high contrast monochrome screen.

## PATIENT DATA CARDS

Built-in card reader provides documentation of patient therapy and treatment outcomes. Five Patient Data Cards come standard with the system. Available card packs of 25.

## THERAPY SYSTEM CART

The optional cart integrates with the unit base to minimize cord and cable hassles.

Six concealed storage bins conveniently house your clinical essentials.

## ORDERING INFORMATION

### System with Cart

- 2793K** – 2 Channel Electrotherapy
- 2797K** – 4 Channel Electrotherapy
- 2791K** – 2 Channel Combination
- 2795K** – 4 Channel Combination

### System

- 2763** – 2 Channel Electrotherapy
- 2786K** – 4 Channel Electrotherapy
- 2760K** – 2 Channel Combination
- 2788K** – 4 Channel Combination

### Accessories

- 2781** – Channel 3/4 Electrotherapy Module
- 2767** – Battery Module
- 2780** – Therapy Cart
- 2790** – Patient Data Management System
- 27465** – Patient Data Cards (25/pack)
- 27381** – 1 cm<sup>2</sup> Applicator
- 27382** – 2 cm<sup>2</sup> Applicator
- 27383** – 5 cm<sup>2</sup> Applicator
- 27384** – 10 cm<sup>2</sup> Applicator
- 27508** – Operator Remote Channels 1 & 2
- 27079** – Operator Remote Channels 3 & 4
- 27468** – Carry Bag

## OPERATOR REMOTES

Conforms comfortably to the hand and allows treatment start, stop or pause at any time as well as change amplitude. (OPTIONAL)

## PATIENT DATA MANAGEMENT (PDM)



Unique optional software program, provides an organized way to store, display, compare and print information captured on Patient Data Cards. (OPTIONAL)

## LEAD WIRE CONNECTIONS

Newly designed lead wires have a rip stop feature and color-coded connectors to easily distinguish between channels.

Patents Pending



ISO 13485 CERTIFIED

4717 Adams Road  
P.O. Box 489  
Hixson, TN 37343 USA  
**1-423-870-2281**  
**1-800-592-7329** USA  
**1-800-361-6661** CANADA  
[www.chattgroup.com](http://www.chattgroup.com)